A decorative graphic on the left side of the slide consisting of two overlapping parallelograms. The front one is blue and the back one is a light greenish-blue. They are positioned diagonally, with the blue one partially covering the green one.

MATLAB to Python Translation: Electron Microscope Tools

COSC 499: Software Engineering



Introduction


- Client: Misa Hayashida, a researcher for the National Research Council
- Project overview: Converting 3 tools from MATLAB to Python
 - qEELS
 - Nanomi Optics
 - Alignment software



UNIVERSITY OF
ALBERTA


qEELS: Use Cases

Name:	qEELS
Description:	Calculate calibrated energy loss from a spectrogram image.
Flow of Events:	<ol style="list-style-type: none">1. The user opens the software.2. The user uploads the spectrogram image.3. The user indicates the location of features on the spectrogram image.4. The user indicates that the software will detect fitted peaks of surface plasmon and bulk plasmon.5. The user requests the calibrated energy loss axis and transfer axis.
Pre-conditions:	<ul style="list-style-type: none">• The user must have the software open.• The user must have an appropriate spectrogram image.• The user has a knowledge of spectrogram images.
Post-Conditions:	The energy loss axis and transfer axis is calculated and displayed in an image.





qEELS: Data Flow

Level 0 DFD



qEELS: Sequence diagram



qEELS: UI Mockup

☒ Bulk Plasmon 1

☐ Bulk Plasmon 1

Width

Detect

X: 921

Y: 251

X:

Y:

60

☐

☐ Bulk Plasmon 1

☐ Bulk Plasmon 1

Width

Detect

X:

Y:

X:

Y:

60

☐

☐ Bulk Plasmon 1

☐ Bulk Plasmon 1

Width

Detect

X:

Y:

X:

Y:

60

☐

Average Pixel

micro rad/pixel upper

micro rad/pixel lower

ev/pixel

Open image

Detect


Save Data

Reset




Nanomi Optics Use Cases

Name:	Nanomi Optics
Description:	Optimize electron microscope optics settings.
Flow of Events:	<ol style="list-style-type: none">1. The user opens the software.2. User inputs their desired parameters for optics3. User request optimized lens settings
Pre-conditions:	<ul style="list-style-type: none">• The user must have the software open.• The user has knowledge of electron microscope settings.
Post-Conditions:	Optimized settings are displayed to the user.







Nanomi Optics Data Flow

Level 0 DFD



Nanomi Optics Sequence Diagrams



Nanomi Optics UI

Settings above sample

Mode: nm

Lens settings (nm):

C1 (nm) ON

C2 (nm) ON

C3 (nm) ON

Settings below sample

Image mode:

☐ Diffraction ☒ Image

Auto Setting:

☐ Objective ☐ Intermediate ☐ Projective ☒ None

Lens settings (nm):

Distance ON

Objective ON

Intermediate ON

Projective ON


Lens Ur: C1 = 0.451495 C2 = 0.658143 C3 = 0.540231

Condensor Aperture: 10

Lens focal length: C1 = 67.290 nm C2 = 22.940 nm C3 = 39.880nm


Magnification: 6167.51

EXTREME beam DIAMETER @ sample = -331.51 nm and convergence SEMI angle = 0.08 mrad






Alignment - Use Cases

Name:	Alignment Software
Description	Align nanoparticle images.
Flow of Events:	<ol style="list-style-type: none">1. The user opens the software.2. The user uploads nanoparticle images.3. The user selects the location of the particle on the image.4. The user requests aligned sequence images from software.
Pre-conditions:	<ul style="list-style-type: none">• The user must have the software open.• The user has appropriate nanoparticle images.• The user has knowledge of nanoparticle images.
Post-Conditions:	An aligned sequence of images is optimized and calculated.




Alignment





Alignment - Data Flow


Level 0 DFD



Alignment - Sequence Diagram



Alignment - Main Window UI



Alignment - Adjust Contrast UI

Contrast Adjustment

Data Range

Minimum:

Maximum:

Window

Minimum:

Maximum:

Width:


Center:

Scale Display Range

☒ Match Data Range

☒ Eliminate outliers ☐

Apply



A histogram showing the distribution of data. The x-axis is labeled from 0 to 1 in increments of 0.1. The y-axis represents frequency. The histogram bars are blue. A red vertical line is drawn at 0.5. A red shaded area covers the entire range from 0 to 1. There are red circular markers at the ends of the x-axis (0 and 1).

Done

Alignment - Automatic Particle Detection and Tracking UI

Automatic Particle Detection and Tracking

Particle detection

	X1	Y1	X2	Y2	IM1	IM2		
1	309	197	389	277	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
2	463	23	543	103	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
3	577	244	657	324	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
4	418	484	498	564	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
5	651	768	731	848	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
6	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
7	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
8	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
9	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
10	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
11	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
12	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track
13	0	0	0	0	1	61	<input type="checkbox"/> Lost	<input type="checkbox"/> Track

Track Selected Particles

Particle color
☒ Black ☐ White

Shift Search Areas

Marker radius (pixel) 20 6

Search area width (pixel) 80 7

Search area height (pixel) 80 7


Up
Left 5 8 Right
Down

RS 9
All RS 10

Pixel info: (X, Y) Intensity 11 1 / 61

Done

Alignment - Manual Tracking Adjustment UI



Alignment - Optimization UI

Optimization

Tilt Angle Specifications

☒ CSV File Start angle:

☐ Constant Step Step angle:

Image Set

[Change Image Set](#)

☒ Fixed rotation angle amount:

☐ Adjust rotation angle amount:


☐ Adjust magnification and rotation:

[Re-perform alignment on images](#)


Accuracy pixels:

/

[Done](#)



Architecture



- Each project will be separate launchable piece of software
- Each project will have a distinct GUI and ENGINE
- The GUI will make calls to the ENGINE
- The ENGINE will not make calls to the GUI
- Common functionality between the three softwares will be factored out into shared Common package



Tech Stack

Programming Languages:

- Python

Python Libraries:

- Tkinter (GUI)
- Matplotlib (Chart Rendering)
- Scipy (Optimization and other algorithms)
- Numpy (Efficient numeric arrays)
- Pillow (Image processing and transforms)
- Flake8 (Linting)
- Pytest (Unit Testing)



Tech Stack

Supporting Technologies:


- Github
- Github Actions (CI/CD)
- Github Projects (Project Management)



Testing Plan


- Set up Github Actions workflows to run all tests automatically
- Pushing to the main branch directly will be blocked
- Merging pull requests without all tests passing first will be blocked
- All pull-requests will require an approving review
- Pull-requests should not be approved unless they include tests which demonstrate the functionality they implement
- All code will be linted with Flake8 automatically

Testing Plan: System Architecture




Testing Plan: Optimization Problems

- Nanomi Optics and the Alignment Software both use Optimization Algorithms
- Optimization output can vary from the original software
- Optimization output may vary run-to-run or machine-to-machine
- Instead of testing for an exact result, use some threshold tests to test the underlying conditions that the software is optimizing for



Testing Plan: Floating-Point Problems

- All three of the softwares use floating-numbers extensively
- Not all numbers can be represented exactly by a floating-point number
- Some loss in precision when operations are performed
- Equality depends on order of operations
- Once again, use thresholds





Testing Plan: Large Output Problems

- Alignment Software performs many operations on large images which much be tested, such contrast adjustment and kernel convolution
- Could be to have an exact copy of the expected output, but that is difficult to hand-craft
- Use smaller sample images for test-cases which will make crafting test-cases easier

